


eUnity Mobile App Searching

To find studies or patients using the eUnity search function, enter patient names or exam-related criteria in any desired combination using the available search fields. Although very similar in use, there are some differences between searching on a web browser versus the mobile app. This document focuses solely on using the search function in the eUnity Mobile App.

Steps for Mobile App Searching

1. Once you've logged into the eUnity app the search window will open with blank fields.




The screenshot shows the eUnity mobile app search interface. At the top, the eUnity logo is visible. Below it, there are several search fields: ACCESSION NUMBER, PATIENT NAME, PATIENT ID, and REFERRING PHYSICIAN, each with a corresponding text input field. Below these are two dropdown menus: DATE RANGE (set to 'Any') and MODALITY (set to 'All'). At the bottom, there is a PATIENT BIRTH DATE section with three input fields for YYYY, MM, and DD. At the very bottom of the screen, there are two buttons: 'Search' and 'Reset'. A small footer at the bottom left reads 'Build #5 0 1 453' and the bottom right reads 'eUnity™ is not certified for diagnostic interpretation on this device.'

Diagnostic Radiology Consultants

2. Enter any combination of patient-related information to help return your desired result.


- a. **Date Range** allows you to filter for exams based on study date.



The screenshot shows the eUnity mobile application interface. At the top, the eUnity logo is visible. Below it, there are several search criteria with input fields: 'ACCESSION NUMBER', 'PATIENT NAME', 'REFERRAL', 'DATE RANGE', 'MODALITY', and 'PATIENT BIRTH DATE'. The 'DATE RANGE' dropdown menu is currently open, displaying a list of options: 'Any' (highlighted in blue), 'Today', 'Yesterday', 'Last 7 Days', and 'Last 30 Days'. Below the search criteria, there are two buttons: 'Search' and 'Reset'. At the bottom of the screen, there is a small disclaimer: 'Build #5.0.1.453 eUnity™ is not certified for diagnostic interpretation on this device.'

Diagnostic Radiology Consultants

b. **Patient Name** allows you to search for your patient by name. When searching by this method, use the format LAST NAME, FIRST NAME. You may also perform partial searches using last name only or wild card searches with first name only using percent signs like so: %name%.




The screenshot shows the eUnity mobile application interface for searching patients. The app has a dark grey background with white text. At the top, the 'eUnity' logo is visible. Below the logo, there are several search criteria, each with a corresponding input field:

- ACCESSION NUMBER**: A text input field.
- PATIENT NAME**: A text input field with the placeholder text 'last name, first name'.
- PATIENT ID**: A text input field.
- REFERRING PHYSICIAN**: A text input field.
- DATE RANGE**: A dropdown menu currently set to 'Any'.
- MODALITY**: A dropdown menu currently set to 'All'.
- PATIENT BIRTH DATE**: Three separate input fields labeled 'YYYY', 'MM', and 'DD'.

At the bottom of the search area, there are two buttons: 'Search' and 'Reset'. Below the buttons, there is a small line of text: 'Build #1.0.1.453 eUnity™ is not certified for diagnostic interpretation on this device.'

Diagnostic Radiology Consultants

c. **Patient ID** allows you to search for your patient by their Medical Record Number (MRN). If you do not know the patient's full MRN but you remember a partial consecutive number you can perform a wild card search with the number between percent signs to return a list of patients with MRNs that match your partial search



The screenshot shows the eUnity mobile application interface for searching patients. The app has a dark grey background with white text. At the top, the eUnity logo is visible. Below the logo, there are several search criteria fields:


- ACCESSION NUMBER: An empty text input field.
- PATIENT NAME: An empty text input field.
- PATIENT ID: A text input field containing the number "123456789".
- REFERRING PHYSICIAN: An empty text input field.
- DATE RANGE: A dropdown menu currently set to "Any".
- MODALITY: A dropdown menu currently set to "All".
- PATIENT BIRTH DATE: Three separate input fields labeled "YYYY", "MM", and "DD", all of which are empty.

At the bottom of the form, there are two buttons: "Search" and "Reset".

Build #5.0.1.453 eUnity™ is not certified for diagnostic interpretation on this device.

Diagnostic Radiology Consultants


d. **Accession Number** allows you to search for your patient by using the exam's Accession Number (ACC#). Partial searches can also be performed by placing numbers between percent signs to return a list of exams matching the wild card search.



The screenshot shows the eUnity mobile application interface for searching. At the top, the eUnity logo and a menu icon are visible. Below are several search criteria fields: 'ACCESSION NUMBER' with a text input containing '%1234%', 'PATIENT NAME', 'PATIENT ID', 'REFERRING PHYSICIAN', 'DATE RANGE' with a dropdown menu set to 'Any', 'MODALITY' with a dropdown menu set to 'All', and 'PATIENT BIRTH DATE' with three input fields for 'YYYY', 'MM', and 'DD'. At the bottom, there are 'Search' and 'Reset' buttons. A small disclaimer at the very bottom reads: 'Suite #5.0.1.453 eUnity™ is not certified for diagnostic interpretation on this device.'

Diagnostic Radiology Consultants


e. **Modality** allows you to search for a specific modality type. By default “ALL” modalities will be returned. Click the drop down and select a modality type to narrow your search results.



The screenshot shows the eUnity mobile application interface. At the top, the eUnity logo is visible. Below it, there are several search fields: 'ACCESSION NUMBER', 'PATIENT', 'PATIENT', 'REFERR', 'DATE R...', and 'MODAL...'. The 'MODAL...' field has a dropdown menu open, showing the following options: 'All' (highlighted in blue), 'CR', 'CT', 'DX', 'MR', and 'MG'. Below the search fields, there are 'Search' and 'Reset' buttons. At the bottom, there is a small disclaimer: 'Suite #5.0.1.453 eUnity™ is not certified for diagnostic interpretation on this device.'

Diagnostic Radiology Consultants

f. Referring Physician lets you search for exams by a specific ordering physician. You can search by just the referring's LAST NAME or the LAST NAME, FIRST NAME. For a quicker results, it helps to also put in a date range or modality.



The screenshot shows the 'Results' screen of the eUnity mobile application. The interface is dark-themed and contains several search criteria fields: 'ACCESSION NUMBER', 'PATIENT NAME', 'PATIENT ID', 'REFERRING PHYSICIAN' (with a placeholder 'LastName, FirstName'), 'DATE RANGE' (set to 'Last 7 Days'), 'MODALITY' (set to 'All'), and 'PATIENT BIRTH DATE' (with fields for 'YYYY', 'MM', and 'DD'). At the bottom, there are 'Search' and 'Reset' buttons. A small disclaimer at the very bottom reads: 'eUnity™ is not certified for diagnostic interpretation on this device.'

3. Click **Search** when you've entered all of your criteria. If you would like to start a new search, click **<Search** at the top left of the search results screen.